

Liver Analysis Synapse 3D

Liver Analysis (CT) is a useful tool for HBP surgeons for pre-operative planning.

The software allows extraction of the liver and nearby vessel regions from contrast-enhanced CT images and displays the results as 3D- mapped images.

Liver Analysis


HPB Surgeons, laparoscopic or open surgery
Useful tool for pre-operative planning

Use as much as you want !

Process, simulate and share, with no limits

Fujifilm unique technology

Liver, Artery, Portal Veins, Hepatic veins, Tumors
Volume calculations
Extraction of vascular territories in liver
Fusion of NM or MRI



Liver Analysis

Benefits

- On site solution
- Identify anatomical variations
- Define your surgery strategy
- Lower surgery time
- Better outcomes and more conservative surgery

