

Dental MPR


Easy Operation to Create Panoramic Image of Alveolar Bone Possible to print DICOM Film, Simulate of implantation

Recommended Image ▶ CT image


Analysis image conditions ▶ Alveolar bone in FOV.

Possible to create Panoramic images along dentition and alveolar bone cross section image. Implant simulation also possible. DICOM film print is corresponded.


Easy Operation



■ Panoramic Image




■ Cross section image of Alveolar Bone



Just click on Alveolar bone then you can have panoramic image and cross section image.

Virtual Implantation



■ Implant Simulation



Can virtually implantation on the image. And can check margin between virtual implant and Mandibular canal

Custom Layout/Print out of Film

■ Screen layout and print layout is editable



Can edit layout as you like.